

# Data Sheet -

## Twist Cap Lancet



### Product Codes:

AMTL-2800

AMTL-3000

Latex Free ✓

Phthalates Free ✓

DEHP Free ✓

### Device Specifications

Device dimensions	Height	31mm
	Width	6mm
	Length	6mm
Weight		0.46g
Material	Housing	Polyethylene
	Lancet	Stainless steel
Lancet	AMTL-2800	28G upto 2.3mm puncture depth
	AMTL-3000	30G upto 2.3mm puncture depth
Sterilisation method		Irradiation by electron beam
Intended use		Capillary blood sampling
Activation		Via lancing device (universal fit)
Puncture site		Finger tips (see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	<b>TBC</b> gC02e

### Packaging Specifications

Pcs/pack		200
Pack weight		100g
Pack dimensions	Height	90mm
	Width	50mm
	Length	70mm
Material		Cardboard (please recycle)
Packs per case		50 (5x inner cases of 10 packs)
Pcs/case		10000
Case weight		6.2kg
Case dimensions	Height	175mm
	Width	280mm
	Length	520mm
Material		Cardboard (please recycle)