

## Data Sheet -

# BV Insulin Pen Needle

### Product Codes:

AMTL-3204

AMTL-3105

Latex Free ✓

Phthalates Free ✓

DEHP Free ✓



## Device Specifications

<b>Device dimensions</b> (includes cover)	Height	30mm
	Width	13mm
	Length	13mm
<b>Weight</b> (includes cover)		1.03g
<b>Material</b>	Protective caps/Needle hub	Polypropylene
	Sealing film	Medical dialysis paper
	Needle	Medical grade stainless steel
<b>Pen needle</b>	AMTL-3204	32G 4mm needle length
		AMTL-3105
	<b>Sterilisation method</b>	
<b>Intended use</b>		Insulin injection
<b>Connection</b>		Standard universal screw thread
<b>Injection site</b>		(see user instructions)
<b>Product carbon footprint</b>	Functional unit: <i>per pack</i>	TBC gC02e

## Packaging Specifications

<b>Pcs/pack</b>		100
<b>Pack weight</b>		116g
<b>Pack dimensions</b>	Height	100mm
	Width	60mm
	Length	105mm
<b>Material</b>		Cardboard (please recycle)
<b>Packs per case</b>		50 (5x inner cases of 10 packs)
<b>Pcs/case</b>		5000
<b>Case weight</b>		7.6kg
<b>Case dimensions</b>	Height	327mm
	Width	230mm
	Length	570mm
<b>Material</b>		Cardboard (please recycle)