

Data Sheet - Band Tourniquet

Product Codes:

AMTL-2803

- Latex Free ✓
- Phthalates Free ✓
- DEHP Free ✓



Device Specifications

Device dimensions	Length	440mm
	Width	25mm
	Thickness	0.5mm
Weight		5.56g
Material	Tourniquet band	TPE
Sterilisation method		Non-sterile
Intended use		To distend veins for phlebotomy
Dispensed		Perforated on a roll
Use site		(see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	TBC gCO2e

Packaging Specifications

Pcs/pack		25
Pack weight		160g
Pack dimensions	Height	110mm
	Width	30mm
	Length	110mm
Material		Cardboard (please recycle)
Packs per case		40
Pcs/case		1000
Case weight		6.5kg
Case dimensions	Height	165mm
	Width	240mm
	Length	460mm
Material		Cardboard (please recycle)