

Data Sheet -

Tri-Plus Adjustable Depth Safety Lancet



Product Codes:

AMTL-2303A

AMTL-2803A

Latex Free ✓

Phthalates Free ✓

DEHP Free ✓

Device Specifications

Device dimensions	Height	53mm
	Width	21mm
	Length	12mm
Weight		2.78g
Material	Housing/Knob/Protective Cap	Polyethylene
	Lancet	Stainless steel
	Spring	Stainless steel & Carbon steel
Lancet	AMTL-2303A	23G 1.3/1.8/2.3mm puncture depth
	AMTL-2803A	28G 1.3/1.8/2.3mm puncture depth
Sterilisation method		Irradiation by electron beam
Intended use		Capillary blood sampling
Activation		Button
Puncture site		Finger tips (see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	TBC gCO2e

Packaging Specifications

Pcs/pack		200
Pack weight		593g
Pack dimensions	Height	110mm
	Width	128mm
	Length	140mm
Material		Cardboard (please recycle)
Packs per case		18
Pcs/case		3600
Case weight		11.1kg
Case dimensions	Height	349mm
	Width	292mm
	Length	398mm
Material		Cardboard (please recycle)