

# Data Sheet - Paediatric Button Tourniquet

## Product Codes:

AMTL-2802

- Latex Free ✓
- Phthalates Free ✓
- DEHP Free ✓



## Device Specifications

<b>Device dimensions</b>	Length	250mm
	Width	22mm
	Thickness	0.7mm
<b>Weight</b>		3.85g
<b>Material</b>	Tourniquet band	TPE
	Tourniquet button	ABS
<b>Sterilisation method</b>		Non-sterile
<b>Intended use</b>		To distend veins for phlebotomy
<b>Dispensed</b>		Singular from pack
<b>Use site</b>		(see user instructions)
<b>Product carbon footprint</b>	Functional unit: <i>per pack</i>	TBC gC02e

## Packaging Specifications

<b>Pcs/pack</b>		100
<b>Pack weight</b>		426g
<b>Pack dimensions</b>	Height	60mm
	Width	140mm
	Length	210mm
<b>Material</b>		Cardboard (please recycle)
<b>Packs per case</b>		10
<b>Pcs/case</b>		1000
<b>Case weight</b>		4.3kg
<b>Case dimensions</b>	Height	225mm
	Width	300mm
	Length	320mm
<b>Material</b>		Cardboard (please recycle)