

Data Sheet -

Pressure Activated Safety Lancet

Product Codes:

AMTL-2618

AMTL-2818

Latex Free ✓

Phthalates Free ✓

DEHP Free ✓



Device Specifications

Device dimensions	Height	40mm
	Width	11mm
	Length	11mm
Weight		1.96g
Material	Housing/Slider/Cap/Cover/Needle holder	Polypropylene
	Needle	Stainless steel
	Spring	Stainless steel & Carbon steel
Lancet	AMTL-2618	26G 1.8mm puncture depth
	AMTL-2818	28G 1.8mm puncture depth
Sterilisation method		Irradiation by electron beam
Intended use		Capillary blood sampling
Activation		Contact with pressure
Puncture site		Finger tips (see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	TBC gCO2e

Packaging Specifications

Pcs/pack		100
Pack weight		210g
Pack dimensions	Height	105mm
	Width	60mm
	Length	110mm
Material		Cardboard (please recycle)
Packs per case		50 (5x inner cases of 10 packs)
Pcs/case		5000
Case weight		12.5kg
Case dimensions	Height	327mm
	Width	230mm
	Length	570mm
Material		Cardboard (please recycle)