

Data Sheet -

BV

Safety Lancet



Product Codes:

AMTL-BV1826 AMTL-BV1830
AMTL-BV1828

Latex Free ✓
Phthalates Free ✓
DEHP Free ✓

Device Specifications

Device dimensions	Height	40mm
	Width	10mm
	Length	7.5mm
Weight		1.10g
Material	Cap/Housing/Slider/Needle holder	Polypropylene
	Needle	Stainless steel
	Spring	Carbon steel
Lancet	AMTL-BV1826	26G 1.8mm puncture depth
	AMTL-BV1828	28G 1.8mm puncture depth
	AMTL-BV1830	30G 1.8mm puncture depth
Sterilisation method		Irradiation by electron beam
Intended use		Capillary blood sampling
Activation		Contact with pressure
Puncture site		Finger tips (see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	TBC gCO2e

Packaging Specifications

Pcs/pack		100
Pack weight		130g
Pack dimensions	Height	65mm
	Width	85mm
	Length	85mm
Material		Cardboard (please recycle)
Packs per case		50
Pcs/case		5000
Case weight		7.1kg
Case dimensions	Height	345mm
	Width	182mm
	Length	435mm
Material		Cardboard (please recycle)