

Data Sheet - Adult Button Tourniquet

Product Codes:

AMTL-2801

- Latex Free ✓
- Phthalates Free ✓
- DEHP Free ✓



Device Specifications

Device dimensions	Length	366mm
	Width	25mm
	Thickness	0.7mm
Weight		5.72g
Material	Tourniquet band	TPE
	Tourniquet button	ABS
Sterilisation method		Non-sterile
Intended use		To distend veins for phlebotomy
Dispensed		Singular from pack
Use site		(see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	TBC gC02e

Packaging Specifications

Pcs/pack		100
Pack weight		664g
Pack dimensions	Height	60mm
	Width	140mm
	Length	240mm
Material		Cardboard (please recycle)
Packs per case		10
Pcs/case		1000
Case weight		7.0kg
Case dimensions	Height	255mm
	Width	300mm
	Length	320mm
Material		Cardboard (please recycle)